

## HISTONE H4 RABBIT PAB

货号: N225449

产品全名: Histone H4 兔多抗

基因符号 Histone H4

**UNIPROT ID:** P62805

**背景:** The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a member of the histone H4 family.

**抗原:** The antiserum was produced against synthesized peptide derived from human Histone H4. AA range:6-55

**经过测试的应用:** WB,IHC-F,IHC-P,ICC/IF,ELISA

**推荐稀释比:** WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 ELISA: 1/10000

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 11 kDa; Observed MW: 11 kDa

**亚型:** IgG


**纯化:** Affinity Purified

**种属反应性:** Human, Mouse and Rat,Monkey


**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Epigenetics and Nuclear Signaling

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Histone H4 in MCF-7 lysates using Histone H4 antibody.



Western blot analysis of Histone H4 in COS7 and HeLa lysates using Histone H4 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---