

HJURP RABBIT MAB

货号: N263366

产品全名: HJURP 兔单克隆抗体

基因符号 FAKTS; URLC9; hFLEG1

UNIPROT ID: Q8NCD3

背景: Centromeric protein that plays a central role in the incorporation and maintenance of histone H3-like variant CENPA at centromeres.

抗原: Recombinant protein of human HJURP

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-4B7

分子量: Calculated MW: 84 kDa; Observed MW: 84 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: -

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using HJURP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of HJURP in HeLa, CHO-K1 lysates using HJURP antibody.