

HLTF RABBIT MAB

货号: N262340

产品全名: HLTF 兔单克隆抗体

基因符号 HIP116; HIP116A; HLTF 1; Hltf; HLTF1; p113; RNF80; SMARC A3; SMARCA 3; SMARCA3; SNF2L3; ZBU1

UNIPROT ID: Q14527

背景: Has both helicase and E3 ubiquitin ligase activities. Possesses intrinsic ATP-dependent nucleosome-remodeling activity; This activity may be required for transcriptional activation or repression of specific target promoters (By similarity).

抗原: A synthetic peptide of human HLTF

经过测试的应用: WB, ICC/IF

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R03-1E4

分子量: Calculated MW: 114 kDa; Observed MW: 114 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of HLTF (green) in HeLa using HLTF antibody, and DAPI (blue)



Western blot analysis of HLTF in HeLa lysates using HLTF antibody.