

## HMGCR RABBIT MAB

货号: N262344

产品全名: HMGCR 兔单克隆抗体

基因符号 LDLCQ3

**UNIPROT ID:** P04035

**背景:** Transmembrane glycoprotein that is the rate-limiting enzyme in cholesterol biosynthesis as well as in the biosynthesis of nonsterol isoprenoids that are essential for normal cell function including ubiquinone and geranylgeranyl proteins.

**抗原:** A synthetic peptide of human HMGCR

**经过测试的应用:** WB,IP

**推荐稀释比:** WB: 1/500-1/1000 IP: 1/20

**种属反应性:** Rabbit

**克隆性:** Rabbit Monoclonal

**克隆编号:** R04-2J3

**分子量:** Calculated MW: 97 kDa; Observed MW: 97 kDa

**亚型:** IgG


**纯化:** Affinity Purified

**种属反应性:** Human,Mouse

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Cardiovascular

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of HMGCR in K562, 3T3, HeLa lysates using HMGCR antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---