

## HNF 4 ALPHA RABBIT MAB

货号: N262346

产品全名: HNF 4 alpha 兔单克隆抗体

基因符号 Hnf4; HNF-4; MODY1; Nr2a1; Tcf14; TCF-14; Hnf4alpha

**UNIPROT ID:** P49698

**背景:** Transcriptionally controlled transcription factor. Binds to DNA sites required for the transcription of alpha 1-antitrypsin, apolipoprotein CIII, transthyretin genes and HNF1-alpha. May be essential for development of the liver, kidney and intestine. Miscellaneous Binds fatty acids.

**抗原:** Recombinant protein of mouse HNF-4-alpha

**经过测试的应用:** WB, IP, ChIP

**推荐稀释比:** WB: 1/500-1/1000 IP: 1/20 ChIP: 1/20

**种属反应性:** Rabbit

**克隆性:** Rabbit Monoclonal

**克隆编号:** R01-7D7

**分子量:** Calculated MW: 53 kDa; Observed MW: 53 kDa

**亚型:** IgG


**纯化:** Affinity Purified

**种属反应性:** Human, Mouse and Rat


**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Cardiovascular

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of HNF4alpha in K562, C6, HeLa lysates using HNF4alpha antibody.



Western blot analysis of HNF4alpha in mouse liver lysates using HNF 4 alpha antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---