

HNF 4 ALPHA RABBIT MAB

货号: N263371

产品全名: HNF 4 alpha 兔单克隆抗体

基因符号 TCF; HNF4; MODY; FRTS4; MODY1; NR2A1; TCF14; HNF4a7;
HNF4a8; HNF4a9; NR2A21; HNF4alpha

UNIPROT ID: P41235

背景: HNF-4-alpha is a transcription factor that belongs to the steroid hormone receptor superfamily. Binds to DNA sites required for the transcription of alpha 1-antitrypsin, apolipoprotein CIII, transthyretin genes and HNF1-alpha.

抗原: Recombinant protein of human HNF-4-alpha

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-4F0

分子量: Calculated MW: 53 kDa; Observed MW: 53 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of HNF4alpha in A549, MCF-7 lysates using HNF 4 alpha antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
