

HNF 4 ALPHA RABBIT MAB

货号: N262345

产品全名: HNF 4 alpha 兔单克隆抗体

基因符号 Hnf4; HNF-4; MODY1; Nr2a1; Tcf14; TCF-14; Hnf4alpha

UNIPROT ID: P49698

背景: Transcriptionally controlled transcription factor. Binds to DNA sites required for the transcription of alpha 1-antitrypsin, apolipoprotein CIII, transthyretin genes and HNF1-alpha. May be essential for development of the liver, kidney and intestine. Miscellaneous Binds fatty acids.

抗原: Recombinant protein of mouse HNF-4-alpha

经过测试的应用: WB, IHC-F, IHC-P, ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-7J6

分子量: Calculated MW: 53 kDa; Observed MW: 53 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Mouse

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of HNF4alpha in Raw264.7 lysates using HNF 4 alpha antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
