

HNRNPA1L2 RABBIT MAB

货号: N262322

产品全名: HNRNPA1L2 兔单克隆抗体

基因符号 hnRNP core protein A1-like 2

UNIPROT ID: Q32P51

背景: Involved in the packaging of pre-mRNA into hnRNP particles, transport of poly(A) mRNA from the nucleus to the cytoplasm and may modulate splice site selection.

抗原: A synthetic peptide of human Heterogeneous Nuclear Ribonucleoprotein (A1-like)

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R05-8G1

分子量: Calculated MW: 34 kDa; Observed MW: 34 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Heterogeneous Nuclear Ribonucleoprotein (A1like) in K562, C6, 3T3, HeLa lysates using HNRNPA1L2 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
