

HOXB4 RABBIT MAB

货号: N262356

产品全名: HOXB4 兔单克隆抗体

基因符号 HOXB4; HOX2F; Homeobox protein Hox-B4; Homeobox protein Hox-2.6; Homeobox protein Hox-2F

UNIPROT ID: P17483

背景: Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

抗原: A synthetic peptide of human HOXB4

经过测试的应用: WB,ICC/IF,IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R05-9A7

分子量: Calculated MW: 28 kDa; Observed MW: 34 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of HOXB4 in K562 lysates using HOXB4 antibody.



Immunocytochemistry analysis of HOXB4 (green) in K562 using HOXB4 antibody, and DAPI (blue).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
