

HRPT2 RABBIT PAB

货号: N225226

产品全名: HRPT2 兔多抗

基因符号 Parafibromin (Cell division cycle protein 73 homolog;
Hyperparathyroidism 2 protein)

UNIPROT ID: Q6P1J9

背景: Tumor suppressor probably involved in transcriptional and post-transcriptional control pathways.

抗原: Synthesized peptide derived from human Parafibromin. at AA range: 51-100

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 61 kDa; Observed MW: 61 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of HRPT2 in mouse liver, mouse brain lysates using HRPT2 antibody.