

HSC70 (5C2) MOUSE MAB

货号: N261372

产品全名: Hsc70 (5C2) 小鼠单抗

基因符号 Heat shock cognate 71 kDa protein (Heat shock 70 kDa protein 8)

UNIPROT ID: P11142

背景: Acts as a repressor of transcriptional activation. Inhibits the transcriptional coactivator activity of CITED1 on Smad-mediated transcription. Chaperone. Component of the PRP19-CDC5L complex that forms an integral part of the spliceosome and is required for activating pre-mRNA splicing.

抗原: Purified recombinant protein expressed in E.coli.

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Mouse

克隆性: Mouse Monoclonal

克隆编号: 5C2-6G4-9A4

分子量: Calculated MW: 71 kDa; Observed MW: 71 kDa

亚型: IgG1


纯化: Affinity Purified

种属反应性: Human,Rat,Mouse


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Hsc70 (5C2) in HeLa, A549, HL-60, U2OS, C6 lysates using Hsc70 (5C2) antibody



Immunohistochemistry analysis of paraffin-embedded Human Breast Carcinoma using Hsc70 (5C2) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.