

HSD17B1 RABBIT MAB

货号: N262364

产品全名: HSD17B1 兔单克隆抗体

基因符号 hydroxysteroid (17-beta) dehydrogenase 1; HSD17; EDHB17; EDH17B2; SDR28C1

UNIPROT ID: P14061

背景: Favors the reduction of estrogens and androgens. Also has 20-alpha-HSD activity. Uses preferentially NADH.

抗原: A synthetic peptide of human HSD17B1

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R09-5K9

分子量: Calculated MW: 35 kDa; Observed MW: 35 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of HSD17B1 in Human placenta lysates using HSD17B1 antibody.