

HSP75 RABBIT PAB

货号: N225742

产品全名: Hsp75 兔多抗

基因符号 HSP75; HSP90L; Trap1

UNIPROT ID: Q12931

背景: TNF receptor-associated protein 1 (TRAP1), also known as HSP75, is a mitochondrial chaperone and ATPase that was originally identified as a protein that interacts with the TNF receptor.

抗原: A synthesized peptide derived from human Hsp75

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 80 kDa; Observed MW: 75 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Hsp75 in K562 lysates using Hsp75 antibody.