

## ICAD RABBIT MAB

货号: N263375

产品全名: ICAD 兔单克隆抗体

基因符号 DFFA; DFF1; DFF45; H13; DNA fragmentation factor subunit alpha; DNA fragmentation factor 45 kDa subunit; DFF-45; Inhibitor of CAD; ICAD

**UNIPROT ID:** O00273

背景: Inhibitor of the caspase-activated DNase (DFF40).

抗原: A synthetic peptide of human ICAD

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-7E7

分子量: Calculated MW: 37 kDa; Observed MW: 45 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ICAD in A549, HL-60, U2OS lysates using ICAD antibody.