

ICAM1 (4F5) MOUSE MAB

货号: N261408

产品全名: ICAM1 (4F5) 小鼠单抗

基因符号 ICAM1; Intercellular adhesion molecule 1; ICAM-1; Major group rhinovirus receptor; CD antigen CD54

UNIPROT ID: P05362

背景: ICAM1 a type I membrane protein of the immunoglobulin superfamily. Is a ligand for the leukocyte adhesion LFA-1 protein (Integrin alpha-L/beta-2) and a Rhinovirus receptor. Typically expressed on endothelial cells and cells of the immune system.

抗原: Purified recombinant human CD54(ICAM-1) protein fragments expressed in E.coli.

经过测试的应用: WB,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200

种属反应性: Mouse

克隆性: Mouse Monoclonal

克隆编号: 4F5-F9-F12

分子量: Calculated MW: 58 kDa; Observed MW: 96 kDa

亚型: IgG2b


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of ICAM1 (4F5) in SW480 using CD54(ICAM1) antibody.



Western blot analysis of CD54(ICAM1) in SW480 and Raji lysates using CD54(ICAM1) antibody.