

ICAM3 RABBIT MAB

货号: N263377

产品全名: ICAM3 兔单克隆抗体

基因符号 CD50; CDW50; ICAM-R

UNIPROT ID: P32942

背景: ICAM proteins are ligands for the leukocyte adhesion protein LFA-1 (integrin alpha-L/beta-2).

抗原: A synthetic peptide of human ICAM3

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-0D3

分子量: Calculated MW: 60 kDa; Observed MW: 110-160 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using ICAM3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of ICAM3 in HeLa, A549, HL-60, U2OS, C6 lysates using ICAM3 antibody.