

IGF-1R(TEPROTUMUMAB) MAB

货号: 28577

产品全名: IGF-1R(Teprotumumab) 单抗

基因符号 CD221, IGF1R, IGF1, JTK13

描述: IGF-1R(teprotumumab) 单抗

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释比: Flow Cyt 1:100

种属反应性: Homo sapiens

亚型: IgG1(E239D, M241L)

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 IGF-1R

成分: Lyophilized from sterile PBS, pH 7.4. 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. IGF-1R protein is highly expressed on the surface of Expi293 cell membrane. Flow cytometry analysis with 1 μ g/mL Anti-IGF-1R (teprotumumab biosimilar) mAb (28577) on Expi293 cells transfected with human IGF-1R (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram), and Isotype antibody on Expi293 transfected with irrelevant protein (Orange histogram)