

I κ B EPSILON RABBIT MAB

货号: N263381

产品全名: I κ B epsilon 兔单克隆抗体

基因符号 NFKBIE; IKBE; NF-kappa-B inhibitor epsilon; NF-kappa-BIE; I-kappa-B-epsilon; I κ B-E; I κ B-epsilon; I κ Bepsilon

UNIPROT ID: O00221

背景: Inhibits NF-kappa-B by complexing with and trapping it in the cytoplasm. Inhibits DNA-binding of NF-kappa-B p50-p65 and p50-c-Rel complexes.

抗原: A synthetic peptide of human I κ B epsilon

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R05-1B2

分子量: Calculated MW: 53 kDa; Observed MW: 48 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of I κ B epsilon in HeLa, Ramos lysates using I κ B epsilon antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
