

IKK EPSILON RABBIT PAB

货号: N226023

产品全名: IKK epsilon 兔多抗

基因符号 IKBKE; IKKE; IKKI; KIAA0151; Inhibitor of nuclear factor kappa-B kinase subunit epsilon; I-kappa-B kinase epsilon; IKK-E; IKK-epsilon; Ikbke; Inducible I kappa-B kinase; IKK-i

UNIPROT ID: Q14164

背景: Phosphorylates inhibitors of NF-kappa-B thus leading to the dissociation of the inhibitor/NF-kappa-B complex and ultimately the degradation of the inhibitor. May play a special role in the immune response. Protects cells against DNA damage-induced cell death.

抗原: A synthetic peptide of human IKKi/IKKe

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF,FC

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 80 kDa; Observed MW: 80 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of IKKi/IKKe in Hela, Raji, Ramos, MCF-7 lysates using IKK epsilon antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
