

IKZF3 RABBIT MAB

货号: N263382

产品全名: IKZF3 兔单克隆抗体

基因符号 IKZF3; ZNFN1A3; Zinc finger protein Aiolos; Ikaros family zinc finger protein 3

UNIPROT ID: Q9UKT9

背景: Transcription factor that plays an important role in the regulation of lymphocyte differentiation. Plays an essential role in regulation of B-cell differentiation, proliferation and maturation to an effector state. Involved in regulating BCL2 expression and controlling apoptosis in T-cells in an IL2-dependent manner.

抗原: A synthetic peptide of human IKZF3

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-0C5

分子量: Calculated MW: 58 kDa; Observed MW: 58 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using IKZF3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of IKZF3 in HL-60, U251 lysates using IKZF3 antibody.