

IL-17B RABBIT PAB

货号: N225221

产品全名: IL-17B 兔多抗

基因符号 IL-17B; interleukin 17B; NIRF; IL-20; ZCYTO7

UNIPROT ID: Q9UHF5

背景: Stimulates the release of tumor necrosis factor alpha and IL-1-beta from the monocytic cell line THP-1.

抗原: Synthetic peptide from human protein at AA range: 1-50

经过测试的应用: IHC-P,ELISA

推荐稀释比: IHC: 1/50-1/100 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: -

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human small intestine using IL-17B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.