

IL-26 RABBIT MAB

货号: N263384

产品全名: IL-26 兔单克隆抗体

基因符号 IL-26; AK155

UNIPROT ID: Q9NPH9

背景: May play a role in local mechanisms of mucosal immunity and seems to have a pro-inflammatory function.

抗原: A synthetic peptide of human IL-26

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R07-3C0

分子量: Calculated MW: 20 kDa; Observed MW: 20 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using IL-26 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of IL-26 in Hela, A549, HL-60, U251 lysates using IL-26 antibody.