

IMPORTIN 9 RABBIT MAB

货号: N261574

产品全名: Importin 9 兔单克隆抗体

基因符号 importin 9; Imp9; Ran-binding protein 9; RanBP9

UNIPROT ID: Q96P70

背景: Nuclear transport receptor that mediates nuclear import of proteins, such as histones, proteasome and actin.

抗原: A synthetic peptide of human Importin 9

经过测试的应用: WB, ICC/IF, IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R03-8D9

分子量: Calculated MW: 116 kDa; Observed MW: 116 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Importin 9/RANBP9 in Jurkat, HeLa lysates using Importin 9/RANBP9 antibody.



Immunocytochemistry analysis of Importin 9 (green) in A549 using Importin 9 antibody, and DAPI (blue).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
