

## INPP4B RABBIT PAB

货号: N225707

产品全名: INPP4B 兔多抗

基因符号 4-bisphosphate 4-phosphatase; Inositol polyphosphate 4 phosphatase type II 105kDa; Inositol polyphosphate 4-phosphatase type II; Inpp4b; MGC132014; Type II inositol 3; 4 bisphosphate 4 phosphatase; Type II inositol-3

**UNIPROT ID:** O15327

**背景:** Catalyzes the hydrolysis of the 4-position phosphate of phosphatidylinositol 3,4-bisphosphate, inositol 1,3,4-trisphosphate and inositol 1,4-bisphosphate.

**抗原:** A synthesized peptide derived from human INPP4B

**经过测试的应用:** WB,IHC-P,FC

**推荐稀释比:** WB: 1/500-1/1000 IHC: 1/50-1/100 FC: 1/50-1/100

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 105 kDa; Observed MW: 105 kDa

**亚型:** IgG


**纯化:** Affinity Chromatography

**种属反应性:** Human, Mouse and Rat

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Signal Transduction

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of INPP4B in MCF-7 lysates using INPP4B antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---