

INTERFERON EPSILON RABBIT PAB

货号: N225178

产品全名: Interferon epsilon 兔多抗

基因符号 IFNE; IFNE1; Interferon epsilon; IFN-epsilon; Interferon epsilon-1

UNIPROT ID: Q86WN2

背景: Type I interferon required for maintaining basal levels of IFN-regulated genes, including 2'-5'-oligoadenylate synthetase, IRF7 and ISG15, in the female reproductive tract.

抗原: The antiserum was produced against synthesized peptide derived from the C-terminal region of human IFNE. AA range:159-208

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 24 kDa; Observed MW: 24 kDa

亚型: IgG

纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Interferon epsilon in Jurkat lysates using Interferon epsilon antibody.