

IQGAP2 RABBIT MAB

货号: N263388

产品全名: IQGAP2 兔单克隆抗体

基因符号 Ras GTPase-activating-like protein IQGAP2

UNIPROT ID: Q13576

背景: Binds to activated CDC42 and RAC1 but does not seem to stimulate their GTPase activity. Associates with calmodulin.

抗原: Recombinant protein of human IQGAP2

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-2C6

分子量: Calculated MW: 181 kDa; Observed MW: 181 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using IQGAP2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of IQGAP2 in A549, HL-60 lysates using IQGAP2 antibody.