

ISLET 1 RABBIT MAB

货号: N262426

产品全名: Islet 1 兔单克隆抗体

基因符号 ISL1; Insulin gene enhancer protein ISL-1; Islet-1

UNIPROT ID: P61371

背景: ISL1 (ISL1 transcription factor, LIM/homeodomain) is a member of the LIM/homeodomain family of transcription factors. It binds to the enhancer region of the insulin gene, among others, and may play an important role in regulating insulin gene expression. It is central to the development of pancreatic cell lineages and may also be required for motor neuron generation. Islet-1 expression defines cardiac progenitor cell populations and is required for normal cardiac development and asymmetry.

抗原: Recombinant protein of human Islet 1

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R05-2C3

分子量: Calculated MW: 39 kDa; Observed MW: 39 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Islet 1 in K562, rat Brain lysates using Islet 1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
