

JUNB RABBIT MAB

货号: N262434

产品全名: JunB 兔单克隆抗体

基因符号 Transcription factor jun-B

UNIPROT ID: P17275

背景: JunB is a basic region, leucine zipper (bZIP) transcription factor belonging to the Jun family that includes c-Jun and JunD. Jun family members homodimerize or heterodimerize with Fos and ATF proteins to form a functional transcription factor AP-1 (activator protein, whose activity is regulated by a variety of physiological and pathological stimuli such as growth factors, infections, and stress signals.

抗原: A synthetic peptide of human JunB

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-9B7

分子量: Calculated MW: 36 kDa; Observed MW: 42,43 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of JunB in K562, C6, 3T3, HeLa lysates using JunB antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
