

## JUNB RABBIT PAB

货号: N225731

产品全名: JunB 兔多抗

基因符号 Transcription factor jun-B

**UNIPROT ID:** P17275

**背景:** JunB is a basic region, leucine zipper (bZIP) transcription factor belonging to the Jun family that includes c-Jun and JunD. Jun family members homodimerize or heterodimerize with Fos and ATF proteins to form a functional transcription factor AP-1 (activator protein, whose activity is regulated by a variety of physiological and pathological stimuli such as growth factors, infections, and stress signals.

**抗原:** A synthesized peptide derived from human JunB

**经过测试的应用:** WB, ICC/IF, IP

**推荐稀释比:** WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 36 kDa; Observed MW: 43 kDa

**亚型:** IgG


**纯化:** Affinity Chromatography

**种属反应性:** Human, Mouse and Rat

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Epigenetics and Nuclear Signaling

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of JunB in HeLa lysates using JunB antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---