

LEUKOTRIENE A4 HYDROLASE RABBIT MAB

货号: N262467

产品全名: Leukotriene A4 Hydrolase 兔单克隆抗体

基因符号 LTA4; LTA4 hydrolase; LTA4H

UNIPROT ID: P09960

背景: Hydrolyzes an epoxide moiety of leukotriene A4 (LTA-4) to form leukotriene B4 (LTB-4). The enzyme also has some peptidase activity.

抗原: A synthetic peptide of human Leukotriene A4 hydrolase

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-6H4

分子量: Calculated MW: 69 kDa; Observed MW: 69 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Leukotriene A4 hydrolase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Leukotriene A4 hydrolase/LTA4H in K562, C6, 3T3, Hela lysates using Leukotriene A4 Hydrolase antibody.