

LLGL1 RABBIT MAB

货号: N263407

产品全名: LLGL1 兔单克隆抗体

基因符号 DLG4; HUGL; LLGL; Lgll; MglI; HUGL1; HUGL-1

UNIPROT ID: Q15334

背景: Cortical cytoskeleton protein found in a complex involved in maintaining cell polarity and epithelial integrity.

抗原: Recombinant protein of human LLGL1

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R05-0B8

分子量: Calculated MW: 115 kDa; Observed MW: 115 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of LLGL1 in A549, U2OS lysates using LLGL1 antibody.