

## LTK RABBIT MAB

货号: N263409

产品全名: LTK 兔单克隆抗体

基因符号 TYK1

**UNIPROT ID:** P29376

背景: Receptor with a tyrosine-protein kinase activity.

抗原: A synthetic peptide of human LTK

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-8A4

分子量: Calculated MW: 92 kDa; Observed MW: 92 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Rat


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: -

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using LTK antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of LTK in A549, HL-60, U2OS, C6 lysates using LTK antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---