

MAPKAP KINASE 3 RABBIT MAB

货号: N263428

产品全名: MAPKAP Kinase 3 兔单克隆抗体

基因符号 MAPKAPK3; MAP kinase-activated protein kinase 3; MAPK-activated protein kinase 3; MAPKAP kinase 3; MAPKAP-K3; MAPKAPK-3; MK-3; Chromosome 3p kinase; 3pK

UNIPROT ID: Q16644

背景: Stress-activated serine/threonine-protein kinase involved in cytokines production, endocytosis, cell migration, chromatin remodeling and transcriptional regulation.

抗原: A synthetic peptide of human MK-3

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-4B6

分子量: Calculated MW: 43 kDa; Observed MW: 43 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of MK3 in A549, HL-60 lysates using MAPKAP Kinase 3 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
