

MATRIN 3 RABBIT MAB

货号: N261553

产品全名: Matrin 3 兔单克隆抗体

基因符号 KIAA0723; Matr3; Matrin3; MPD2; VCPDM

UNIPROT ID: P43243

背景: May play a role in transcription or may interact with other nuclear matrix proteins to form the internal fibrogranular network. In association with the SFPQ-NONO heteromer may play a role in nuclear retention of defective RNAs.

抗原: A synthetic peptide of human Matrin 3

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R09-5E3

分子量: Calculated MW: 95 kDa; Observed MW: 125 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling


储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human colon cancer tissue using Matrin 3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Matrin 3 (green) in HL-60 using Matrin 3 antibody, and DAPI (blue).



Western blot analysis of Matrin 3 in Jurkat, Hela lysates using Matrin 3 antibody.