

## MDB3 RABBIT PAB

货号: N226051

产品全名: MDB3 兔多抗

基因符号 MBD3; Methyl-CpG-binding domain protein 3; Methyl-CpG-binding protein MBD3

**UNIPROT ID:** O95983

背景: Does not bind DNA by itself. Recruits histone deacetylases and DNA methyltransferases. Acts as transcriptional repressor and plays a role in gene silencing.

抗原: A synthetic peptide of human MBD3

经过测试的应用: WB, ICC/IF, FC, IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 33 kDa; Observed MW: 33 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of MDB3 in A375, A431, mouse brain lysates using MDB3 antibody.



Western blot analysis of MDB3 in C6 lysates using MDB3 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---