

## MKNK2 RABBIT PAB

货号: N225456

产品全名: MKNK2 兔多抗

基因符号 MKNK2; GPRK7; MNK2; MAP kinase-interacting serine/threonine-protein kinase 2; MAP kinase signal-integrating kinase 2; MAPK signal-integrating kinase 2; Mnk2

**UNIPROT ID:** Q9HBH9

背景: Serine/threonine-protein kinase that phosphorylates SFPQ/PSF, HNRNPA1 and EIF4E.

抗原: Synthesized peptide derived from the Internal region of human Mnk2.

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 52 kDa; Observed MW: 52 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: -

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of MKNK2 in HeLa lysates using MKNK2 antibody.