

## MMAE(8B4) RABBIT MAB

货号: 28570

产品全名: MMAE(8B4) 兔单克隆抗体

基因符号 N.A.

描述: MMAE antibody(8B4) 兔单克隆抗体

背景: N.A.

经过测试的应用: ELISA

推荐稀释比: ELISA 1/5000-10000

种属反应性: Rabbit

亚型: Rabbit IgG

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: N.A. MMAE

成分: Lyophilized from sterile PBS, pH 7.4. 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at  $-80^{\circ}\text{C}$  (Avoid repeated freezing and thawing).

ELISA assay to evaluate Anti-MMAE antibody  
0.2 $\mu\text{g}$  Human IgG-MMAE per well




Figure1. Elisa plates were pre-coated with IgG-MMAE (0.2 $\mu\text{g}$ /per well). Serial diluted anti-MMAE monoclonal antibody (28570) solutions were added, washed, and incubated with secondary antibody before Elisa reading. From above data, the EC50 for anti-MMAE monoclonal antibody binding with IgG-MMAE is 2.274ng/ml.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---