

MMAE(8C4) RABBIT MAB

货号: 28571

产品全名: MMAE(8C4) 兔单克隆抗体

基因符号 N.A.

描述: MMAE antibody(8C4) 兔单克隆抗体

背景: N.A.

经过测试的应用: ELISA

推荐稀释比: ELISA 1/5000-10000

种属反应性: Rabbit

亚型: Rabbit IgG

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: N.A. MMAE

成分: Lyophilized from sterile PBS, pH 7.4. 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

ELISA assay to evaluate Anti-MMAE antibody
0.2 μg Human IgG-MMAE per well




Figure1. Elisa plates were pre-coated with IgG-MMAE (0.2 μg /per well). Serial diluted anti-MMAE monoclonal antibody (28571) solutions were added, washed, and incubated with secondary antibody before Elisa reading. From above data, the EC50 for anti-MMAE monoclonal antibody binding with IgG-MMAE is 1.356ng/ml.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
