

MMP9 RABBIT MAB

货号: N262530

产品全名: MMP9 兔单克隆抗体

基因符号 92 kDa gelatinase; 92 kDa type IV collagenase; Gelatinase B; GELB

UNIPROT ID: P41245

背景: Matrix metalloproteinase that plays an essential role in local proteolysis of the extracellular matrix and in leukocyte migration (By similarity).

抗原: Recombinant protein of mouse MMP9

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-3A5

分子量: Calculated MW: 81 kDa; Observed MW: 92 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of MMP9 in rat lung lysates using MMP9 antibody.