

## MOG RABBIT PAB

货号: N225879

产品全名: MOG 兔多抗

基因符号: BTN6; BTNL11; MOGIG2; NRCLP7

**UNIPROT ID:** Q16653

**背景:** Mediates homophilic cell-cell adhesion (By similarity). Minor component of the myelin sheath. May be involved in completion and/or maintenance of the myelin sheath and in cell-cell communication.

**抗原:** A synthesized peptide derived from human Myelin oligodendrocyte glycoprotein

**经过测试的应用:** WB, IHC-P, ICC/IF

**推荐稀释比:** WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 28 kDa; Observed MW: 28 kDa

**亚型:** IgG

**纯化:** Affinity Chromatography

**种属反应性:** Human, Mouse and Rat

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Neuroscience

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Myelin oligodendrocyte glycoprotein in mouse brain lysates using MOG antibody.