

MSP RABBIT MAB

货号: N261585

产品全名: MSP 兔单克隆抗体

基因符号 serine threonine kinase 3/4; KRS1; MST2/KRS2; MST1; YSK3; TIIAC

UNIPROT ID: Q13043

背景: Stress-activated, pro-apoptotic kinase which, following caspase-cleavage, enters the nucleus and induces chromatin condensation followed by internucleosomal DNA fragmentation.

抗原: A synthetic peptide of human Serine/threonine-protein kinase 4

经过测试的应用: WB, ICC/IF, IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R09-3C1

分子量: Calculated MW: 56 kDa; Observed MW: 60 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of MSP/MST1 in K562, rat Brain, C6, 3T3, HeLa lysates using MSP/MST1 antibody.



Immunocytochemistry analysis of MSP (green) in K562 using MSP antibody, and DAPI (blue).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
