

## MST4 RABBIT MAB

货号: N262542

产品全名: MST4 兔单克隆抗体

基因符号 MASK; MST4

**UNIPROT ID:** Q9P289

背景: Mediator of cell growth (PubMed:11641781, PubMed:17360971). Modulates apoptosis (PubMed:11641781, PubMed:17360971). In association with STK24 negatively regulates Golgi reorientation in polarized cell migration upon RHO activation (PubMed:27807006).

抗原: A synthetic peptide of human MST4

经过测试的应用: WB,ICC/IF,IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-2H9

分子量: Calculated MW: 47 kDa; Observed MW: 52 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Rat


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of MST4 (green) in HeLa using MST4 antibody, and DAPI (blue)



Western blot analysis of MST4 in K562, C6 lysates using MST4 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---