

MTAP RABBIT MAB

货号: N263437

产品全名: MTAP 兔单克隆抗体

基因符号 BDMF; MSAP; DMSFH; LGMBF; DMSMFH; c86fus; HEL-249

UNIPROT ID: Q13126

背景: This gene encodes an enzyme that plays a major role in polyamine metabolism and is important for the salvage of both adenine and methionine. The encoded enzyme is deficient in many cancers because this gene and the tumor suppressor p16 gene are co-deleted. Multiple alternatively spliced transcript variants have been described for this gene, but their full-length natures remain unknown. [provided by RefSeq, Jul 2008]

抗原: Recombinant protein of human MTAP

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R07-418

分子量: Calculated MW: 31 kDa; Observed MW: 31 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of MTAP in A549, MCF-7, U251, HL-60 lysates using MTAP antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
