

MTTFA RABBIT PAB

货号: N225918

产品全名: mtTFA 兔多抗

基因符号 MtTF1; mtTFA; TCF6; TCF6L1; TCF6L2; TCF6L3; TFAM;
Transcription factor 6

UNIPROT ID: Q00059

背景: Binds to the mitochondrial light strand promoter and functions in mitochondrial transcription regulation. Required for accurate and efficient promoter recognition by the mitochondrial RNA polymerase.

抗原: A synthesized peptide derived from human mtTFA

经过测试的应用: WB,IHC-P,ICC/IF,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 29 kDa; Observed MW: 24 kDa

亚型: IgG

纯化: Affinity Chromatography

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Tags & Cell Markers

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of mtTFA in MCF-7 lysates using mtTFA antibody.