

MUSCARINIC ACETYLCHOLINE RECEPTOR M2 RABBIT MAB

货号: N262553

产品全名: Muscarinic Acetylcholine Receptor M2 兔单克隆抗体

基因符号 CHRM2; HM2; AChR; Acm2

UNIPROT ID: P08172

背景: The muscarinic acetylcholine receptor mediates various cellular responses, including inhibition of adenylate cyclase, breakdown of phosphoinositides and modulation of potassium channels through the action of G proteins. Primary transducing effect is adenylate cyclase inhibition.

抗原: A synthetic peptide of human Muscarinic Acetylcholine Receptor 2

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-9B9

分子量: Calculated MW: 52 kDa; Observed MW: 52 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience&Muscarinic Acetylcholine Receptor

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Muscarinic Acetylcholine Receptor 2/CM2 in K562, C6, HeLa lysates using Muscarinic Acetylcholine Receptor M2 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
