

MVK RABBIT MAB

货号: N263442

产品全名: MVK 兔单克隆抗体

基因符号 MVK; Mevalonate kinase; MK

UNIPROT ID: Q03426

背景: Catalyzes the phosphorylation of mevalonate to mevalonate 5-phosphate, a key step in isoprenoid and cholesterol biosynthesis.

抗原: A synthetic peptide of human MVK

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-6D3

分子量: Calculated MW: 42 kDa; Observed MW: 42 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: -

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of MVK in HeLa, A549, HL-60, U251, U87-MG lysates using MVK antibody.