

## MYC MAB

产品名称: MYC 小鼠单克

货号: 28001

规格: 100  $\mu$ L

**Symbol:**C-Myc

描述:MYC 小鼠单抗

背景: The Anti-MYC antibody recognizes any protein containing Myc Tag.

**Immunogen:** CEQKLISEEDL

经过测试的应用: ELISA, WB, IF, IHC

推荐稀释比: ELISA 1:2000-1:5000WB 1:1000-1:5000IF 1:100-1:200

**Concentration:** 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG


纯化: Purified from ascites

**Preservative:** No

成分: PBS (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Recognizes protein contain Myc .

储存条件: Store at  $-20^{\circ}C$ . Avoid repeated freezing and thawing



蛋白质印迹: 蛋白质印迹: 重组的c-Myc标签融合蛋白, 用c-Myc鼠单克隆抗体(1:1000)稀释), 蛋白上样量为50ng进行蛋白印迹分析