

MYH1 (9C7) MOUSE MAB

货号: N261219

产品全名: MYH1 (9C7) 小鼠单抗

基因符号 fast skeletal muscle; MYH3; HEMHC; Myosin-1; Myosin heavy chain 2x; MyHC-2x; Myosin heavy chain IIx/d; MyHC-IIx/d; Myosin heavy chain; skeletal muscle; adult 1

UNIPROT ID: P12882

背景: Muscle contraction.

抗原: Synthetic Peptide of Myosin Heavy Chain

经过测试的应用: ICC/IF

推荐稀释比: IF: 1/50-1/200

种属反应性: Mouse

克隆性: Mouse Monoclonal

克隆编号: 9C7-4A6-5H10

分子量: -

亚型: IgG1


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat, Fruitfly, Nematode

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunofluorescence analysis of MYH1 (9C7) in Immunofluorescence Staining of nematode tissue using MYH1 (9C7) antibody.