

MYH10 RABBIT MAB

货号: N263458

产品全名: MYH10 兔单克隆抗体

基因符号 NMMHCB; NMMHC-IIB

UNIPROT ID: P35580

背景: Cellular myosin that appears to play a role in cytokinesis, cell shape, and specialized functions such as secretion and capping.

抗原: A synthetic peptide of human non-muscle Myosin IIB/MYH10

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R05-6D8

分子量: Calculated MW: 229 kDa; Observed MW: 229 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of non-muscle Myosin IIB/MYH10 in A549, HL-60, U251, U87-MG lysates using non-muscle Myosin IIB/MYH10 antibody.



Western blot analysis of non-muscle Myosin IIB/MYH10 in 3T3 lysates using MYH10 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
