

NAV1.7 RABBIT PAB

货号: N225151

产品全名: Nav1.7 兔多抗

基因符号 SCN9A; NENA; Sodium channel protein type 9 subunit alpha; Neuroendocrine sodium channel; hNE-Na; Peripheral sodium channel 1; PNI; Sodium channel protein type IX subunit alpha; Voltage-gated sodium channel subunit alpha Nav1.7

UNIPROT ID: Q15858

背景: This gene encodes a voltage-gated sodium channel which plays a significant role in nociception signaling. Mutations in this gene have been associated with primary erythralgia, channelopathy-associated insensitivity to pain, and paroxysmal extreme pain disorder.

抗原: The antiserum was produced against synthesized peptide derived from human SCN9A. AA range:651-700

经过测试的应用: WB,IHC-P,ELISA

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 226 kDa; Observed MW: 220 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Nav1.7 in various lysates using Nav1.7 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
